


A meandering introduction into the world of light microscopy

Josh Lawrimore, PhD

8/12/2025

The trilemma: The challenges of using light



Magnification: Refraction of light using lenses


- We, particularly Japan and Germany, have gotten really good at making lenses
 - Chromatic aberration corrections
 - Spherical aberration corrections
 - Anti-reflective coatings
- Objectives typically go up to magnification of 100x
- Not much point in going past 100x when using light due to the diffraction limit

Diffraction Limit

⚠ Equations on this page! ⚠

$$\text{Resolution} = \frac{0.61\lambda}{NA}$$

$$NA = n \sin \theta$$




- n = refractive index of media (air 1.0, water 1.33, oil 1.515)
- θ = half-angle of the maximum cone of light that the lens can capture
- λ = Wavelength of the light
- NA = Numerical Aperture of objective

<https://www.microscopyu.com/microscopy-basics/numerical-aperture>

<https://www.microscopyu.com/microscopy-basics/refractive-index-of-refraction>


<https://www.microscopyu.com/microscopy-basics/resolution>

The trilemma: The challenges of using light




A brief introduction to contrast


- Classic microscopy techniques typically use stains to provide contrast
- For example:
Hematoxylin and
Eosin (H&E)



Fluorescence: A New Kind of Contrast



Diffraction Pattern




This center peak looks like a Gaussian distribution

Breaking the diffraction limit: 2014 Nobel Prize in Chemistry


- Fluorophores are point sources of light that will create a diffraction pattern
- If you could **isolate single fluorophores in time**, you could fit the diffraction pattern with a 2D Gaussian

Stimulated Emission Depletion Microscopy



- Creates a sub-diffraction-limited **excitation zone**
- Repress the surrounding emission by creating a **depletion zone**
- Scan the region to reconstruct a super resolution image

PALM: Photoactivated Localization Microscopy



- Use a special PA-GFP that uses UV light to convert inert fluorophores to an active one
- Only use small amount of UV light to convert a small sample
- Intentionally bleach active PA-GFPs
- Then start new activation cycle

Cost-Benefit of Enhanced Resolution

